# Fact Sheet SMA-ED

## NEW STYRENE-MALEIC ANHYDRIDE (SMA) POLYMERS

Detergent-free system for structural and functional studies of membrane proteins

### What is it?

- Polymer-forming planar lipid bilayer nanodiscs for membrane protein reconstitution
- Direct, detergent-free reconstitution of membranes
- Stable in a broad pH range and in proximity to divalent cations
- Good for structural and functional studies of membrane proteins
- Easy nanodisc size control adjust by changing polymer:lipid ratio


#### Why use it?

• Part of Anatrace's new SMA portfolio


anatrace

CALIBRE SCIENTIFIC COMPANY


- Stable at all pH conditions (excluding those between 5 and 7) [1]
- Found to be tolerant to Ca<sup>2+</sup> and Mg<sup>2+</sup> for all tests at pH 3.5 (from 10 to 200 mM) [1]



pH stability of SMA-ED, adapted from [2]. Stabilizing pH range shown in green.



"And the other states



DLS profiles showing the different sized nanodiscs obtained by varying the lipid-to-polymer weight ratio for SMA-ED [1].

#### Background

- Styrene:Maleic anhydride 1:1
- Molecular Weight ~7.8 kDa
- Solubility (Water) ≤20%
- pH (1% in water) 2.0

Item #

SMA-ED 1 G

SMA-ED 500 MG

SMA-ED 250 MG

Structure of SMA-ED, with ethylene diamine functionalization

**Ordering Information** 

SMA-ED

SMA-ED

SMA-ED

Supporting Documentation: SDS • CoA


Description

UOFM

ΕA

ΕA

FA



**UOM 2021** 

\$560

\$310

\$175

#### **Applications**

- Studies of membrane proteins in native lipid environments
  - Structural studies of large membrane proteins and complexes
- Studies of protein:lipid interactions
- Work with proteins unstable in detergents
- Solubilization of membrane proteins at low or high pH conditions

#### References

 Ravula, T. et al. Effect of polymer charge on functional reconstitution of membrane proteins in polymer nanodiscs. Chem Commun 54, 9615–9618 (2018).
Ravula, T., Hardin, Nathaniel. Z., Mauro, G. M. D. & Ramamoorthy, A. Styrene maleic acid derivates to enhance the applications of bio-inspired polymer based lipid-nanodiscs. Eur Polym J 108, 597–602 (2018).

[3] Ravula, T., Hardin, N. Z. & Ramamoorthy, A. Polymer nanodiscs: Advantages and limitations. Chem Phys Lipids 219, 45–49 (2019).

[4] Chen, A., Majdinasab, E. J., Fiori, M. C., Liang, H. & Altenberg, G. A. Polymer-Encased Nanodiscs and Polymer Nanodiscs: New Platforms for Membrane Protein Research and Applications. Frontiers Bioeng Biotechnology 8, 598450 (2020).

Phone : 800.252.1280 Fax : 419.740.6630



customerservice@anatrace.com www.anatrace.com

